

Compact Orbital Shaker for CO₂ incubator

SANYO

MIR-S100C

**Orbital shaker
for CO₂ incubator**

**Designed for
high humidity use**

Revolution speed

40 to 200 rpm



MIR-S100C
Installed in SANYO
CO₂ incubator

Compact orbital shaker for culture incubation in high humidity conditions

- Usable in a CO₂ incubator.
- Controller and shaker sections are separated to prevent humidity from affecting electrical components.
- Copper alloy stainless steel is used for the shaker section.
- Revolution speed is from 40 to 200 rpm.
- Best suited to operation in a SANYO CO₂ incubator.



Copper-enriched
Stainless Steel Platform

Incubator, culture vessels, clamps and platform are optional.

MIR-S100C

Compact orbital shaker for culture incubation in high humidity. Easy installation in SANYO CO₂ incubator.

Features

• **Specially designed for incubation even in high-humidity conditions**

Controller section and shaker section are independently constructed so electrical components are not affected by humidity. Cable connection socket is waterproof.

• **Installs easily in following SANYO CO₂ incubators.***

MCO-20AIC/19AIC/19AIC(UV)/19M/19M(UV)/18AIC/18AIC(UV)/18AC
(shaking section: W409 x D275mm)

Functions

Structural design

Antibacterial copper alloy stainless steel is used for the shaker section. Also, the shaker section is equipped with a small window that can be used for spraying disinfectant inside.

Soft start

Gradually accelerates until reaching the preset number of revolutions so that samples inside flasks and test tubes are not spilled or foamed. Three levels of acceleration speed can be selected.

Timer

Indicates time of operation from 0 to 99.9 hours in increments of 0.1 hour (6 minutes).

Auto recovery after power interruption

Automatically recovers after power interruption and indicates that power interruption has occurred.

Speed error alarm

Indicates an error when the number of revolutions does not reach or exceeds the preset value.

Motor overload protection circuit

Detects over-current and automatically stops operation. Indicates an error when overload is applied to the motor.

Maintenance-required indicator

Automatically alerts you when maintenance is required or when bearings and belts need replacing after operational hours exceed 10,000.

Options

Platform



MIR-38PFT
(uses metric screws)
Threaded hole diameter: 10 – 32UNF
Threaded hole pitch (longitudinal): 27.5 mm (1.08")
Threaded hole pitch (lateral): 27.5 mm (1.08")

MIR-38PFN
(uses US standard screws)
Threaded hole diameter: 10 – 32UNF
Threaded hole pitch (longitudinal): 34.85 mm (1.37")
Threaded hole pitch (lateral): 20.65 mm (0.81")

*The platform also accommodates clamps and accessories manufactured by other companies.


Clamps for flasks

MIR-51FC	for 50ml flask
MIR-101FC	for 100ml flask
MIR-201FC	for 200ml flask
MIR-501FC	for 500ml flask
MIR-1001FC	for 1000ml flask
MIR-2001FC	for 2000ml flask





MIR-301FC
for 300ml flask

Spring rack
MIR-38SP
for various tubes




Test tube rack
MIR-204SC
for various tubes





Control panel



RPM indicator lamp ————

Display of revolution speed and time ————

START/STOP key ————

TIME indicator lamp ————

Numeric value shift key ————

Digit shift key ————

ENT key ————

SET key ————


SELECT key ————

Specifications

Model	MIR-S100C	
External dimensions (W x D x H)	Shaker	409 x 275 x 118 (mm) / 16.1 x 10.8 x 4.6 (inch)
	Controller	365 x 275 x 93 (mm) / 14.4 x 10.8 x 3.7 (inch)
Exterior	Shaker	Body: Copper alloy stainless steel
	Controller	Body: Galvanized steel with baked-on finish, acrylic coating
Optional platform dimensions (W x D)	380 x 250 (mm) / 15.0 x 9.8 (inch)	
Optional platform	Copper alloy stainless steel	
Net weight	16kg (35.3 lbs.)	
Shaking motion	Circling	
Amplitude (Circling diameter)	20mm	
Revolution	40 to 200 rpm	
Allowable load weight	5kg (11 lbs.)	
Flask intake capacity	Shaker	100ml Erlenmeyer flask x 15, 500ml Erlenmeyer flask x 5 (MIR-38PFT)
	Controller	100ml Erlenmeyer flask x 11, 500ml Erlenmeyer flask x 6 (MIR-38PFN)
Display	LED digital 3-digit display (Speed/Time switching)	
Motor	DC brushless motor	
Applied environmental condition	Shaker	In incubator, 0 to 50°C
	Controller	0 to 40°C without dew condensing
Accessories	Shaker	Rubber plug with a cable hole x 2, Connection cable (3m) x 2,
	Controller	Shaker stand x 1, Rubber foot x 4, Screw x 4, Double-deck platform

Caution: SANYO guarantees the product under certain warranty conditions.
SANYO in no way shall be responsible for any loss of content or damage to content.
* SANYO is not responsible for damage to other than SANYO-branded incubators when used with this Shaker.
• Appearance and specifications are subject to change without notice.

Dimensions



SANYO Electric Co., Ltd., Biomedical Division, Gunma Factory is certified for:



Quality management system: ISO9001
Medical devices quality management system: ISO13485

SANYO Electric Co., Ltd., Biomedical Division, Gunma Factory is certified for:

Environmental management system: ISO14001 (JACO)

Distributed by:



SANYO Electric Co., Ltd.
Biomedical Division
<http://biomedical.sanyo.com/>